

enta [®] ceram polychrome		enta [®] ceram		enta [®] ceram solution	
Anteriors Up	Anteriors Low*	Posteriors		Anteriors Up	Anteriors Low*
L42C	M42	E2, E4, 34S, H4		S10	S2
M42C	M42	E2, E4, 34S, H4		S20	S4-S6
M45C	M45	E4, 34S		S30	S4-S6
S45C	M45	E4, 34S		S40	S6
M48C	M48	E5, 34		S50	S6
M51C	M48	E5, 34			
M45CA	M45	E4, 34S			
M45FA	M45	E4, 34S			
S48FA	S48	E4, 34			
L46CS	M42	E2, E4, 34S, H4			
L46FS	M42	E2, E4, 34S, H4			

* When using Optiform, model 30, 34 or 34S (Lingualized Occlusion), use a 2 to 3 mm wider set of lower anteriors.

		enta [®] ceram classic						
		Anteriors			Posteriors			
		Up	Low with E/H	Low with Optiform	20° [E]	10° [H]	New Criterion	Optiform
Oval	110	62-70	70-72	E2	H4	292	34S-30	
	120	62	70-72	E4-E5		312	34	
	130	70	70-72	E2	H4	292	34S-30	
	140	72	76	E4		312	34	
	150	76	78	E2-E4	H4	292 - 312	34S-30	
Rectangular	210	64	72	E4		312	34	
	220	64	72	E4		312	34	
	230	74	76	E5			34	
	240	66-80	80-84	E4		312	34	
	250	76	84	E4-E5		312	34	
	260	82	84-88	E5			34	
	270	82	84-88	E5			34	
	280	86	86-88	E5			34	
Triangular	310	60	62	E2	H4	292	34S	
	320	68	76	E2	H4	292	34S	
	330	64	72	E4		312	34	
	340	74	76	E4-E5		312	34	
	350	70	70-72	E2	H4	292	34S-30	
	360	64-68	72	E2-E4	H4	292 - 312	34S-30	
	370	68-76	76	E4-E5		312	34	
	380	72-76	76-84	E4		312	34	
Quare	410	62	70-72	E4		312	34-34S	
	420	74	76-84	E4		312	34-34S	
	430	64	72	E4		312	34-34S	
	440	80	82	E5			34	
	450	76-82	84	E5			34	
	460	74	76-84	E4		312	34-34S	

enta[®]cryl
polychrome

enta[®]cryl

enta[®]cryl
basic

Anteriors Up	Anteriors Low*	Posteriors	Anteriors Up	Anteriors Low*	Posteriors
M39C	M39	E2, 34S, H4	M39C	M39	E2, 34S, H4
M42C	M42	E2, E4, 34S, H4	M42C	M42	E2, E4, 34S, H4
M45C	M45	E4, 34S	M45C	M45	E4, 34S
M48C	M48	E5, 34	M48C	M48	E5, 34
S48C	S48	E4, H4	S48C	S48	E4, 34
M51C	M48	E5, 34	M51C	M48	E5, 34
M42F	M42	E2, E4, 34S, H4	M42F	M42	E2, E4, 34S, H4
M45F	M45	E4, 34S	M45F	M45	E4, 34S
M48F	M48	E5, 34	M48F	M48	E5, 34

enta[®]cryl
classic

Posteriors Up	Anteriors Low*	Posteriors
V10	L40-L45	P7-ST4
V20	L25-L40	P4-ST2
V30	L25-L40	P4
ST10	L20-L25	P2
ST20	L45-L55	P7-ST4
ST30	L25-L40	P4-ST2
ST40	L20-L25	P2-P4
ST50	L20-L25	P2-P4
O10	L25	P7-ST4
T10	L10-L20	P4-ST2
T20	L25	P4-P2
T30	L10-L20	P2
T40	L20	P2-P4
D5	L10	P2
D15	L40	P4-ST2
D25	L25	P4-ST2
D28	L40	P4
D30	L50	P4
D35	L40-L50	P4
D40	L52-L45	P4
D45	L52-L55	ST4

Cross bite situation: P5 posteriors

* When using Optiform, model 30, 34 or 34S (Lingualized Occlusion), use a 2 to 3 mm wider set of lower anteriors.

contact

enta b.v.
Brouwerijbaan 10
4615 AA Bergen op Zoom
The Netherlands

Phone: +31 (0)164.23.74.54
Fax: +31 (0)164.23.53.12
E-Mail: info@enta.nl
Internet: www.enta.nl

