

# VERGLEICHSTABELLE FÜR PORZELLAN ZAHNFORMEN UND FARBEN / COMPARISON CHART FOR PORCELAIN MOULDS AND SHADES

## Vitapan Lumin Vacuum vs. Candulor PhysioSet CT Porcelain

Formen OK / Upper moulds

O

Vita	Candulor CT
513	201/302
182	332/343
612	352/452
180	302/201
565	303
185	353
742	334/204

△

Vita	Candulor CT
627	302
734	232
622	302
732	302/232
737	263
773	203/264
843	304

□ □

Vita	Candulor CT
632	412
190	423
722	412/413
187	412/413
743	212/424
854	203
843	304

Formen UK / Lower moulds

Vita	Candulor CT
539	53
639	53
760	53/81
740	81/82
730	53/82
55	83
843	304

Farben / Shades

Vita Classic	Candulor CT
A1	J0
A2	J1
A3	J2
A3.5	J3
A4	M4/M5
B1	J0
B2	M1
B3	M3
B4	M4
C1	901
C2	S2
C3	S2
C4	S3
D2	S1
D3	902
D4	903/S2

## Ivoclar Vivoperl PE vs. Candulor PhysioSet CT Porcelain

Formen OK / Upper moulds

O

Ivoclar PE	Candulor CT
A32	452
A54	343/353
A36	343
A56	353/334

△

Ivoclar PE	Candulor CT
A44	202
A41	232
A24	232/233
A22	302
A24B	233/232
A42	202
A25	203/304
A27	203/204

□ □

Ivoclar PE	Candulor CT
A11	412
A13	412
A12	263/412
A66	424
A14	203/413
A15	425/264
A68	264/424
A17	204/224

Formen UK / Lower moulds

Ivoclar PE	Candulor CT
A4	81
A3	81/53
A2	81
A6	82
A5	82
A7	82
A9	83
A8	83

Farben / Shades

Vivoperl PE	Candulor CT
110 / 01	J0
120 / 1A	J1
140 / 1C	J2
230 / 1E	J3

130 / 2A	M1
210 / 2B	M2
320 / 5B	M3
340 / 3E	M4
240 / 2C	M5

410 / 4A	S1
510 / 6D	S2
530 / 3C	S3
D3	902
D4	903/S2



**Candulor AG**  
Pünten 4 8602 Wangen Switzerland  
Tel. +41(0) 44 805 9000 Fax +41(0)44 805 9090  
[www.candulor.com](http://www.candulor.com) [candulor@candulor.ch](mailto:candulor@candulor.ch)

# VERGLEICHSTABELLE FÜR PORZELLAN ZAHNFORMEN UND FARBEN / COMPARISON CHART FOR PORCELAIN MOULDS AND SHADES

## Dentsply Bioblend vs. Candolor PhysioSet CT Porcelain

Formen OK / Upper moulds



Dentsply Bioblend	Candolor CT
62D	332
62G	453
65G	334



Dentsply Bioblend	Candolor CT
42C, 42D	302
43D	202
44E, 42D	262
43F, 42F	203
45F	233
42G	304
45H	204
42H	234



Dentsply Bioblend	Candolor CT
22C, 32B, 52C, 21C	423
11D, 13E, A84, 55D, 21D	412
32E, 22E, 24F, 54F, 55F	212
12E, 75E, 21X, 21E	424
31F, 11G, 22G, 25G	413
32G	425
12G, 22H	213
11H, 21H, 21G, 21J	224

Formen UK / Lower moulds

Dentsply Bioblend	Candolor CT
C, D, E, F, H	80
K	81
M, N	52
P	82
O, R	53
U, S	83
K1	54
Z	84

Farben / Shades

Dentsply Bioblend	Candolor CT
100, 102, 104	J1
108, 110, 112	J2
109, 118	J3
104	M1
114	M3
116	M4
118	M5
112, 106	S1
113	S2

## Enta Ceram vs. Candolor PhysioSet CT Porcelain

Formen OK / Upper moulds



Enta Ceram	Candolor CT
110	332
120	262
130	353/423
140	342/352
150	452/343



Enta Ceram	Candolor CT
210	332
220	232
230	263/303
240	303/343
250	264
260	264
270	233/304
280	413



Enta Ceram	Candolor CT
310	201
320	332
330	262
340	202
350	262/263
360	412
370	203
380	264/333



Enta Ceram	Candolor CT
410	332/262
420	342/412
430	412
440	212
450	212/452
460	413

Formen UK / Lower moulds

Enta Ceram	Candolor CT
60	-
62	-
64	80
66	80/52
68	80
70	52
72	53
74	81/82
76	82
78	82/53
80	82
82	82/83
84	83
86	83
88	84/54

Farben / Shades

Enta	Candolor CT
A1	J0
A2	J1
A3	J2
A3.5	J3
A4	M4/M5
B1	J0
B2	M1
B3	M3
B4	M4
C1	901
C2	S2
C3	S2
C4	S3
D2	S1
D3	902
D4	903/S2

## Enta Solution vs. Candolor PhysioSet CT Porcelain

Formen OK / Upper moulds



Enta Solution	Candolor CT
S10	332/202
S20	334/203
S30	203/413
S40	233
S50	453/304

Formen UK / Lower moulds

Enta Solution	Candolor CT
S2	53/81
S4	82/53
S6	83/54



**Candolor AG**  
Pünten 4 8602 Wangen Switzerland  
Tel. +41(0) 44 805 9000 Fax +41(0)44 805 9090  
www.candolor.com candolor@candolor.ch

Ivoclar and Vivopearl PE are registered trademarks of Ivoclar Vivadent Schaan, Liechtenstein.  
Vita and Lumin Vacuum are registered trademarks of Vita, Bad Sackingen, Germany.  
Dentsply and Bioblend are registered trademarks of Dentsply, York PA.  
Enta Ceram and Enta Solution are registered trademarks of Enta, Bergen op Zoom, Netherlands.