



**Cross Linked Multi Layered  
Acrylic Teeth**  
Cross Linked Multi Layered  
Kunststof Elementen



**Super-C<sup>®</sup>**

# anterior Super-C<sup>®</sup>

tanden



A B C

9.1 7.8 41.3

8.0 5.1 33.1

**13**



9.6 8.0 41.7

8.2 5.2 33.0

**25**



9.7 8.0 43.1

8.6 5.5 35.5

**35**



A B C

9.9 7.9 42.9

8.4 5.4 35.3

**36**



9.7 8.2 43.9

8.6 5.5 35.5

**56**



A B C

10.0 8.3 45.5

9.0 5.7 36.1

**78**



11.3 8.9 47.5

11.7 5.6 37.4

**98**



12.5 9.7 51.4

10.6 6.1 41.4

**99**



A B C

10.0 8.0 42.9

8.6 5.4 35.4

**66**




10.2 8.1 44.8

9.3 5.7 35.9

**76**





Cross Linked Multi Layered  
Acrylic Teeth

**Super-C<sup>®</sup>**

A B C

10.8 8.4 45.4

9.3 5.1 34.4

**57**



10.3 8.5 45.3

9.7 5.6 38.0

**77**



10.5 8.6 46.1

8.6 5.6 36.2

**88**



posteriors **Super-C<sup>®</sup>**

kiezen

A B C

9.5 8.5 46.3

8.6 5.6 35.9

**85**



A 30.4  
B 33.4  
C 8.5  
D 7.35

**13**



A 29.2  
B 31.0  
C 8.8  
D 8.3

**41**



A 32.0  
B 34.0  
C 9.1  
D 8.6

**44**



A 34.5  
B 35.9  
C 9.4  
D 8.4

**45**



A 29.5  
B 31.3  
C 8.1  
D 7.1

**11**




A 30.2  
B 32.0  
C 9.0  
D 7.7

**12**




**Mould 13** model 13

- A 9.1 — Length of central without collar (mm)
- B 7.8 — Width of central (mm)
- C 41.3 — Width of 6 anteriors on flat plane (mm)




**Mould 13** model 13

- A 8.0 — Length of central without collar (mm)
- B 5.1 — Width of central (mm)
- C 33.1 — Width of 6 anteriors on flat plane (mm)




**Mould 11** model 11

- A 29.5 — The average width of a 1x4 set of upper posteriors
- C 8.1 — The average depth of the upper left first molar



**Mould 11** model 11

- B 31.3 — The average width of a 1x4 set of lower posteriors
- D 7.1 — The average depth of the lower left first molar



**combination**  
combinatie

anteriors tanden	posteriors kiezen	
13	41	12
25	41	11
35	41	44
98	45	
99	45	
36	41	11
56	12	
66	12	
76	41	11
77	44	13
88	44	13
78	41	
85	44	13
57	44	13

**shades**  
kleuren

A1	B0	C2	D3
A2	B1		
A3	B2		
A3.5			
A4			

ISO  
22112



SC 358-2



017-1

CE 0373



Volg ons / Follow us:



Facebook Instagram Youtube LinkedIn



enta® b.v.  
Brouwerijbaan 10  
4615 AA Bergen op Zoom  
The Netherlands

Phone +31 (0) 164.237.454  
Fax +31 (0) 164.235.312  
E-mail info@enta.nl  
Internet www.enta.nl