

a bit more variety...

enta[®]cryl



Vormenkaart
Mould chart
Formenkarte
Carte de formes
acrylic



enta[®]

anterior upper		upper anterior teeth 21 moulds			Central L1	Central W1	Set of W2
V10		11,60	9,35	49,30			
V20		10,60	8,50	45,50			
V30		10,60	8,80	46,70			
ST10		10,50	8,05	44,00			
ST20		12,40	10,20	54,50			
ST30		10,90	8,75	46,50			
ST40		9,85	8,60	45,40			
ST50		10,75	8,50	45,00			
Q10		11,50	8,60	46,00			
T10		10,25	8,00	43,50			
T20		10,40	8,45	46,00			

anterior upper		Central L1	Central W1	Set of W2
T30		9,60	8,80	45,50
T40		8,70	7,85	41,50
D5		9,00	7,40	40,00
D15		10,25	8,25	43,50
D25		11,90	9,20	47,00
D28		10,10	8,70	47,50
D30		10,20	8,10	45,30
D35		11,30	8,80	45,60
D40		10,80	8,60	48,20
D45		12,70	9,00	49,20

anterior's upper
upper anterior teeth
9 moulds

Central	Central	Set of
L1	W1	W2



anterior's lower
lower anterior teeth
5 moulds

Central	Central	Set of
L1	W1	W2



anterior lower	lower anterior teeth 10 moulds	Central	Central	Set of
		L1	W1	W2
L10		8,85	4,95	32,00
L20		9,10	5,20	33,50
L25		9,80	5,40	35,50
L40		9,65	5,70	37,20
L45		11,85	5,90	37,50
L48		8,80	4,90	33,00
L50		9,40	5,40	35,50
L51		10,50	5,20	36,30
L52		10,40	5,50	37,30
L55		11,40	5,80	38,20

NL

pritidenta b.v. beschikt over een grote ervaring wat betreft productontwikkeling en de vervaardiging van porselein en acryl prothese-elementen. pritidenta b.v. werkt nauw samen met verschillende universiteiten in binnen- en buitenland.

Enta® tanden en kiezen voldoen aan de hoogste eisen wat betreft sterkte, slijtvastheid, anatomische vormgeving en natuurgetrouwe karakterisering en kleur.

Enta® producten voldoen aan alle internationale standaarden en normen. De A1-D4 benaming correspondeert met de Vita kleurring. Dit is bedoeld als oriëntatie, niet als vergelijking. Vita is geen geregistreerd merk van pritidenta b.v.

GB

pritidenta b.v. has had large experience concerning research and development and the production of porcelain and acrylic teeth. pritidenta b.v. cooperates closely with several dental universities in the Netherlands and abroad.

Enta® teeth meet the highest requirements in respect of strength, grind ability, anatomical forms, natural characterization and shades. Enta® products meet all international standards.

D

pritidenta b.v. verfügt über eine große Erfahrung hinsichtlich Produktentwicklung und Herstellung von Prothesenzähnen aus Kunststoff und Keramik. pritidenta b.v. arbeitet eng mit Universitäten in den Niederlanden und aus dem Ausland zusammen.

Enta® Zähne genügen höchsten Anforderungen im Hinblick auf ihre Stärke, Schleifbarkeit, anatomische Formgebung und naturgetreue Charakterisierung und Farben. Enta® Produkte entsprechen allen internationalen Normen.

F

pritidenta b.v. a une grande expérience en ce qui concerne la recherché, le développement et la production de dents en porcelaine et en résine. pritidenta b.v. collabore avec des universités aux Pays-Bas et à l'étranger.

Les dents Enta® répondent aux critères, les plus exigeants en ce qui concerne la robustesse, la possibilité de rectifier, la forme anatomique et la reproduction des couleurs naturelles. Les produits Enta® satisfont à toute les normes internationales.

posterior lower	upper	OPTIFORM 3 moulds	1/2 Set of		1. Molars	
			W3	W4	W3	W4
30			29,50	8,20		
			30,00	7,80		
34			32,80	9,50		
			34,00	8,50		
34S			32,40	9,20		
			33,00	8,20		

- L1 Length 1:**
Length of the central incisors without the cervical [in mm]
- W1 Width 1:**
Width of the central incisors [in mm]
- W2 Width 2:**
Width of the total set [in mm]
- W3 Width 3:**
Width of the half set [in mm]
- W4 Width 4:**
Width of the first molar [in mm]



ISO 3336: 1993
ISO 4824: 1993

- All moduls are illustrated actual size
- Available in A0 - D4, including B0
- The dimensions are in millimetres (mm)
- The A1 - D4 designations correspond to the Vita shade guide
- This is meant to be a guide, not a match
- Vita is not a registered trademark of pritidenta b.v.

posteriors
lower upper

1/2 Set of
W3 W4
1. Molars



enta[®] cryl

basic & Verilux posteriors



posterior
lower upper

posterior teeth
8 moulds


1/2 Set of
W3 W4
1. Molars
W4



a bit more personal...

Your enta contact

Do you still have any questions or would you like to arrange for a personal consulting session? We look forward to your call or e-mail.



Enta b.v.
Brouwerijbaan 10
4615 AA Bergen op Zoom
The Netherlands

Phone: +31 (0)164.23.74.54
Fax: +31 (0)164.23.53.12
E-mail: info@enta.nl
Internet: www.enta.nl